

Opilex®

OPILEX® WOUND DRESSING™ Roll

OPILEX® WOUND DRESSING Roll is a disposable non-woven adhesive dressing for wound protection. OPILEX® WOUND DRESSING Roll has an easy peeling technology based on a hypoallergenic adhesive, and also is designed to prevent secondary infection at the application regions with dual safety technology.

OPILEX® WOUND DRESSING Roll is consist of high air permeable pad that rapidly absorb exudate fluid, and provides a stable tissue healing environment. It uses a soft non-woven fabric for protection of sensitive skin and air permeability.



How to use WOUND DRESSING




Description

Name	OPILEX® WOUND DRESSING Roll
Classification	Quasi-drug
Purpose of use	Wound protect
Instruction for use	1 film on application region(PRN)
Ingredient	Non-woven adhesive tape/Polyethylene mesh embo pad
Manufacture	OPILEX Co., LTD.
Country of manufacture	Republic of Korea

Standard specification

Product No.	Size	Packing unit	CTN
OPW-0404	4 cm × 4 cm	100 PCS	100 EA



Contact for exclusive distributor

Phone : +82-51-728-8202 / Fax : +82-51-727-8203
E-mail : opilex@opilex.co.kr

OPILEX CO., LTD.

Customer service




Address : 4-8, Gadong 2-gil, Jeonggwang-eup, Gijang-gun,
Busan, Republic of Korea
Phone : +82-1577-5264 / E-mail : opilex@opilex.co.kr
Web : www.opilex.co.kr



OPILEX[®] WOUND DRESSING[™] Roll

Dressing innovative technologies of OPILEX Co., LTD.

Dressing accuracy of treatment
Dressing adhesive does not stick to hand
Dual removal technology of release paper
Secondary infection protective technology



OPILEX CO., LTD.