

Opilex®


OPILEX® VISIBLE DRESSING™

오피렉스 비저블 드레싱은 상처 부위 및 환부 등의 보호 용도로 사용되는 폴리우레탄 점착 드레싱입니다.

비저블 드레싱은 **Dual Safe**(2중 PET이 형지 설계)와 **Easy Peeling**(쉬운 박리)기술을 적용하여 점착제가 손에 닿지 않아 발생할 수 있는 감염(Infection)을 방지하여 보다 위생적입니다. **Air Way**(통기성 엠보 패드)를 적용한 패드는 우수한 삼출물 흡수와 통기 효율을 높여 보다 안정된 치료 환경을 조성합니다. 저자극 점착제를 사용한 비저블 드레싱의 Film은 **Water Block**(방수)기능으로 수분, 먼지, 세균 등을 차단하여 2차 감염을 방지하며, **Ultra Thin**(초박형 두께)로 굴곡진 부위의 적용에도 아주 얇게 밀착됩니다. 투명한 재질의 폴리우레탄 필름은 의료 행위 후 직접적인 모니터링이 가능합니다.



How to use VISIBLE DRESSING




Description

제품명	오피렉스 비저블드레싱
제품구분	의약외품
효능·효과	상처부위, 환부 등의 보호
용법·용량	환부에 1일 1회 내지 수회 1매씩 갈아붙인다.
전성분	폴리우레탄점착필름, 패드(부직포/폴리에틸렌망)
제조원	주식회사 오피렉스/부산광역시 북구 호열로 111 부산지식산업센터 109, 110호
원산지	MADE IN KOREA

Standard specification

Product No.	Size	Packing unit	CTN	EDI금액
OPVD-4555	4.5 cm × 5.5 cm	50 EA	2,000 EA	M3202134
OPVD-6070	6 cm × 7 cm	50 EA	2,000 EA	M3202134
OPVD-0910	9 cm × 10 cm	50 EA	2,000 EA	M3203134
OPVD-0915	9 cm × 15 cm	50 EA	2,000 EA	M3204134
OPVD-0920	9 cm × 20 cm	50 EA	1,500 EA	M3205134
OPVD-0925	9 cm × 25 cm	50 EA	1,000 EA	M3205134
OPVD-0930	9 cm × 30 cm	50 EA	1,000 EA	M3205134
OPVD-0935	9 cm × 35 cm	50 EA	1,000 EA	M3205134



고객문의

주소 : 부산광역시 기장군 정관읍 가동2길 4-8(1F)

Tel : 051-728-8202 / Fax : 051-727-8203

E-mail : opilex@opilex.com / Web : www.opilex.co.kr

Customer service

Address : 4-8, Gadong 2-gil, Jeonggan-eup, Gijang-gun, Busan, Republic of Korea

Phone : +82-1577-5264 / E-mail : opilex@opilex.co.kr / Web : www.opilex.co.kr

Opilex® 오피렉스® Opi® 오피아이® 는 주식회사 오피렉스의 의료용품 브랜드입니다.
Opilex® 오피렉스® Opi® 오피아이® is the medical product brand of Opilex.



OPILEX CO., LTD.



OPILEX[®] VISIBLE DRESSING[™]

Dressing innovative technologies of OPILEX Co., LTD.

Dressing accuracy of treatment
Dressing adhesive does not stick to hand
Dual removal technology of release paper
Secondary infection protective technology



OPILEX CO., LTD.