

Opilex®

OPI® INNOCARE IV™

오피아이 이노케어 I.V는 의료 기구용 클립, 폴리우레탄필름으로 구성되어 있으며, 인체조직을 제외한 카테터 및 수액 튜브 등의 처치 시 일시적으로 고정하거나 잡아주기 위해 사용되는 I.V 드레싱입니다.

저 자극 점착제를 사용한 오피아이 이노케어 I.V는 Dual Safe(2중 PET이 형지 설계)와 Easy Peeling and Safe Removal(쉬운 박리와 안전한 제거)기술을 적용하여 점착제가 손에 닿지 않으며, 환부의 부착 및 박리에 의한 감염 및 피부손상(Infection and Skin Damage)을 방지하여 보다 위생적이고 안전합니다.



How to use INNOCARE IV



Description

제품명	오피아이 이노케어 아이브이
품목명	의료기구용 클립
모델명	INNO-A
제품구분	일회용 의료기기
멸균방법	EO GAS 멸균
제조원	주식회사 오피렉스/경상남도 김해시 소망길 88 메디컬실용화센터 606호
원산지	MADE IN KOREA

Standard specification

Product No.	Size	Packing unit	CTN	EDI급여
OPIV-0693	6cm × 9.3cm	50 PCS	1,600 EA	M6740234



고객문의

주소 : 부산광역시 기장군 정관읍 가동2길 4-8(1F)

Tel : 051-728-8202 / Fax : 051-727-8203

E-mail : opilex@opilex.com / Web : www.opilex.co.kr

Customer service

Address : 4-8, Gadong 2-gil, Jeonggwan-eup, Gijang-gun, Busan, Republic of Korea

Phone : +82-1577-5264 / E-mail : opilex@opilex.co.kr / Web : www.opilex.co.kr

OPILEX CO., LTD.



Opilex® 오피렉스® Opi® 오피아이® 는 주식회사 오피렉스의 의료용품 브랜드입니다.
Opilex® 오피렉스® Opi® 오피아이® is the medical product brand of Opilex.



OPI[®] INNOCARE IV[™]

Dressing innovative technologies of OPILEX Co., LTD.

Dressing accuracy of treatment
Dressing adhesive does not stick to hand
Dual removal technology of release paper
Secondary infection protective technology



OPILEX CO., LTD.